



Masters Company Inc.
Creative Solutions for Water Management

FILTER PAPER



Precise and accurate testing can be aided with filter paper. Available in a variety of styles and thickness.

| CAT. NO. | PRODUCT DESCRIPTION | |
|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|
| 611 | Pre-Treatment Filtration Test Kit. Consists of: 30 ml Plastic Syringe, Filter Holder Assembly, Whatman No. 5 or Equivalent Grade Filter Paper, 2.5 cm, & Sample Pre-Treatment Filtration Instructions (100 Sheets/Box) | Each |
| 760 | Filter Paper, Whatman No. 1 or Equivalent, 11.0 cm (Fast Filter Rate) (100/Box) | Each |
| 761 | Filter Paper, Whatman No. 5 or Equivalent, (For Pressure Filter Assemblies) 2.5 cm (Slow Filter Rate) (100/Box) | Each |
| 762 | Filter Paper, Whatman No. 5 or Equivalent, 11.0 cm (Slow Filter Rate) (100/Box) | Each |
| 763 | Filter Housing, Pressure Filter Syringe Fitment, (For use with 2.5 cm Filter Circles, No. 5 or No. 0.45 um Filtration Papers) | Each |
| 764 | Filter Paper, Whatman No. 5 or Equivalent 12.5 cm (Slow Filter Rate) (100/Box) | Each |
| 764A | Filter Paper, 0.45 micron, 2.5 cm (100/Box) (Ultra Filtration) (For use in Pressure Filter Syringe Units) | Each |
| 769 | Filter Paper, Whatman No. 1 or Equivalent, 12.5 cm (Fast Filter Rate) (100/Box) | Each |